

# Meeting #10

11/11/21

Omar Luna

# Deliverables

- Complete function to create epochbm\_dict
- Document authored functions
- Work on repo readme file

# Methodology and Learnings

- Did more research on vhdr files
  - BrainVision documentation
  - Timestamp placed in different location
    - Unable to pull out timestamp
- Fixed comments and uploaded `pd_epoch_dict.py`
- Worked on readme file

# Results

---

bmarr

Functions:

`pd_epoch_dict(tobii_data, epoch_dict, eeg_ts)`

Function to return a dictionary containing the epochs as keys and a list of boundaries and data as the values

Libraries:

- numpy

ADELE Dependencies:

- decisionTree\_epochDetection
- getEpochbm\_dict
- js\_pupil\_diameter
- read\_eeg

ADELE Dependents:

- none

Inputs:

`tobii_data`

Pandas DataFrame containing the pupil diameter data with timestamp values as index.

`epoch_dict`

Dictionary with predetermined epochs and boundaries.

`eeg_ts`

The initial time of recording of the Tobii Pro Glasses 2 as a string.

Outputs:

`epochbm_dict`

Dictionary containing the integer value of the epochs as keys and a list containing the epoch boundaries and pupil diameter data as the values.

Description:

The `pd_epoch_dict()` function converts the `eeg_ts` timestamp into the `datetime.datetime` format. Then, the function iterates through the `epoch_dict` values to convert the boundary values into the `datetime.datetime` format and appends the values to the `l_epoch_ts` list. For each `epoch_dict` value, the epoch boundaries are used to determine which data values from `tobii_data` within that epoch. Those data values are appended to the `l_list` list. A tuple is created containing the epoch boundaries from `l_epoch_ts` and the data values from `l_list`. The tuple created is appended to the list `l_tps`. Once all of the boundary and data tuples are stored, the `epochbm_dict` dictionary is created by using the `zip()` function to pair the `epoch_dict` keys with the corresponding tuples from `l_tps`. The function then returns `epochbm_dict`.

```
json_to_dataframe(participant_filename, data_filename)
```

Function to parse JSON files for pupil diameter data and start time.

Libraries:

- JSON
- pandas
- numpy
- datetime.datetime

ADELE Dependencies:

- None

ADELE Dependents:

- pd\_epoch\_dict

Inputs:

participant\_filename

JSON file produced by running an experiment with the Tobii Pro Glasses 2 containing the information regarding the individual using the device as well as information about recorded data.

data\_filename

JSON file produced by running an experiment with the Tobii Pro Glasses 2 containing the data values recorded during the experiment.

Outputs:

df\_final

The pandas DataFrame containing the left, right, and average pupil diameters as columns and the timestamps of the data as index values.

Description:

The json\_to\_dataframe() function retrieves the initial time of the recording from the participant\_filename using the get\_ptimestamp() function and the data recorded from the experiment using the get\_data() function. Once the data is stored, it is used to create a dataframe that is reshaped through DataFrame.pivot() which sets the timestamp values as the index and the left and right pupil diameters as the columns. The initial time of recording is then converted into a timestamp value and used to adjust the values in the DataFrame. A final DataFrame with the converted index values and the left, right, and average pupil diameters is returned.

Sub-Functions:

-get\_ptimestamp(filename)

Function that returns the initial time of the recording from the participant\_filename.

Inputs:

filename

JSON file produced by running an experiment with the Tobii Pro Glasses 2 containing the information regarding the individual using the device as well as information about recorded data.

Outputs:

dt\_timestamp

The initial time of recoding as a datetime.datetime object

### 6.3.2 Marker File

BrainVision Data Exchange Marker File, Version 1.0

[Common Infos]

Codepage=UTF-8

DataFile=example.eeg

[Marker Infos]

; Each entry: Mk<Marker number>=<Type>,<Description>,<Position in data points>,  
; <Size in data points>, <Channel number (0 = marker is related to all channels)>  
; Fields are delimited by commas, some fields might be omitted (empty).

; Commas in type or description text are coded as "\".

Mk1=New Segment,,1,1,0,20130909192513857877

---

Brain Vision Data Exchange Marker File, Version 1.0

[Common Infos]

Codepage=UTF-8

DataFile=2020\_06\_04\_T05\_U00T\_EEG01.eeg

[Marker Infos]

Mk1=Stimulus,35842,1,0,

Mk2=Stimulus,35843,2178,0,

Mk3=Stimulus,35842,2424,0,

### Header file

[Comment]

Recording Start Time

2020-06-04 23:45:19.217

Channel	Name	Impedance (kOhm)	Offset (mV)
0	Fp1	16	5
1	Fp2	16	8